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Sonic tip for activating rinsing liquids

1. Indication
 In combination with the tip holder SF1981/SFQ1981 and the clamping nut 587, the sonic tip SF65 (rinsing tip) is used for sound-based activation of rinsing and disinfecting solutions previously inserted into the root canal and for removing bacterial biofilms. The vibration induced in the instrument tip is transferred to the liquid in the root canal, causing the current hereby generated to remove the smear layer. In addition, this enables the rinsing liquid to penetrate even the tiniest of cavities (e.g. dentinal tubules), thus ensuring an effective disinfection of the root canal.

2. Contraindication
 The sonic tips are not suitable for a very narrow and/or extremely curved root canals, or only in the straight sections of such canals. Failure to observe this may lead to an inadvertent tearing/tearing of the root canal or the preparation of a "via falsa".

3. Warnings
 The useful life of the sonic tips is limited. They therefore have to be used with utmost care.
 Attention:
 - Risk of fracture! Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with utmost care as even minor mechanical damage can provoke the fracture of the instruments.
 - Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this (for example by using a rubber dam).
 - Premature wear! Work in narrow or curved root canals can damage the rinsing tip. Likewise, excessive contact pressure can lead to damage to the rinsing tip or tip holder. Inspect instruments after every use for signs of wear or damage. Exchange if necessary.
 - Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced.
 - Excessive pressure in axial direction! Excessive pressure in axial direction can lead to the creation of a "via falsa". Therefore, make sure to insert the sonic tip gently into the canal, no further than the first point of resistance. To control the penetration depth, consult a radiograph and only insert the tip into the straight or slightly curved section of the canal. The tip has been provided with a laser mark for orientation. Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons. Monitor and record carefully how many times each instrument has been used (5–8 canals at most, when used in extremely curved canals, discard after one use, if necessary)
 Safety control:
 Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger. For additional safety, load the tips mechanically with approx. 10 N (1 kg), without function. Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection.

4. Start-up
 The instruments are supplied in a non-sterile condition and have to be prepared prior to first use.
 The assembled sonic tip can be used:
 - In the Komet sonic hand piece SF1LM/S
 - In the SONICflex™ hand piece made by the co. KaVo (Series 2003N/L/X/LX)
 - In scalars made by co. W&H (i.e. Series Proxo® ZA-55/L/LM/M/LS or Proxo® ST ZE-55 RM/BC)
 - In the SIROAIR L provided by co. Dentsply Sirona
 Make sure to turn off the water supply of the dental treatment unit. The tip changer with torque SF1975 is used for inserting and removing the tip holder. Screw the chosen tip holder into the hand piece with the mentioned tip changer, turning the tool in clockwise direction. To remove tip holder, place tip changer on the tip and turn anticlockwise.
 Note: Please observe the directions in the instructions supplied with your sonic hand piece.
 Secondly, to assemble the rinsing tip ref. SF65 screw the clamping nut (Fig. 5). Now screw on the clamping nut (Fig. 5). Now screw the shank end of the rinsing tip or two pitches in a clockwise direction. The recess in the clamping nut must remain open (Fig. 3)! Now insert it up with the axis of the tip holder (Fig. 4). Pull rinsing tip forward a little until it lies against the clamping nut. Hold the preassembled rinsing tip diagonally downwards, to prevent it from falling out of the clamping nut (Fig. 5). Now screw on the clamping nut by turning it in a clockwise direction until the stop position has been reached (Fig. 6). To secure the connection, tighten the clamping nut with the torque wrench or the tip changer. To loosen the clamping nut, turn in anticlockwise direction with the torque wrench or tip changer. Completely unscrew the clamping nut. Slide the rinsing tip towards the back to remove it from the clamping nut. Forward the product components to the reprocessing venue.

5. Proper use
 - First, treat the root canal in the usual manner by means of an activated preparation (e.g. Komet F360, AlphaKite or EasyShape).
 - An existing rinsing protocol can be expanded by adding the sound established system of the root canal. The final rinsing cycle is carried out with NaOCl in the usual concentration. The tip SF65 can be integrated in a rinsing protocol according to the following example:
 a. Rinse the canal with NaOCl and activate for 15 seconds.
 b. Rinse with a new NaOCl solution.
 c. Carry out a further rinsing and activating sequence.
 d. Remove NaOCl with EDTA and activate for 15 seconds.
 e. Remove EDTA with NaOCl and activate for 15 seconds.
 The rinsing solution has to be inserted into each canal to be treated by means of a cannula. To ascertain the insertion depth, determine the length of the root canal with the help of a radiograph or with mechanical or electric measuring methods. Do not allow the tips to jam in the canal.
 - Once the rinsing tip has been inserted up to 2–3 mm short of the apex (attention: for narrow and curved canals, see section 3 "warnings"), actuate it at power setting 1. To facilitate this, the SF65 is provided with helpful depth marks at 16–18 mm and 20–22 mm. The use of power settings 2 and 3 is not permitted when working with these tips. During use, move tip up and down performing gentle stroke movements. This prevents the rinsing tip from concentrating on the same spot for too long.
 - The provision of plenty of cooling or changing the rinsing solution during the treatment will improve the efficiency.
 - After use, dry rinsing tip thoroughly (e.g. with alcohol) and continue treatment in the usual manner.

6. Maintenance
 Place disassembled sonic tip or individual components (see point 4) in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

7. Reprocessing
 The instruments are supplied in a non-sterile condition and have to be prepared prior to first use.
 Manual pre-cleaning (prior to mechanical and manual reprocessing):
 1. Remove tip holder, clamping nut and rinsing tip from cleaning/disinfection tank (Fräsator) immediately before mechanical reprocessing.
 2. Remove stubborn contamination with a nylon brush under running water, turning the tip holder, rinsing tip and clamping nut constantly.
 3. Check patency of the cooling passage of the tip holder by guiding the nozzle cleaner through the passage. Replace instrument if the passage is obstructed.

Mechanical reprocessing (validated):

1. Rinse tip holder, clamping nut and rinsing tip thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfecter.
 2. Unscrew the nozzle on the injector bar of the washer/disinfecter or pull it out. Then screw or push the external thread of the rinse adapter SF1978/SFQ1978 in the silicon holders on the injector bar. Now screw the tip holder in the internal thread at the top of the rinse adapter (Fig. 7). Observe the directions in the instructions for use provided with the SF1978/SFQ1978 (Fig. 3). The rinsing tip and the clamping nut can be processed separately in the thermo disinfecter, inside suitable containers.
 3. Put chemical detergent into the washer/disinfecter, following the indications on the label and the instructions of the manufacturer of the washer/disinfecter.
 4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place following the A₀ value and observing national provisions (prEN/ISO 15883).
 5. On completion of the cycle remove tip holder, clamping nut and rinsing tip from the washer/disinfecter and dry (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the tip holder or clamping nut.
 6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.
 Manual reprocessing (alternative):
 1. Place tip holder and clamping nut onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforation, rinse again with detergent/disinfectant.
 2. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip and the individual components constantly.
 3. Rinse tip holder, clamping nut and rinsing tip thoroughly with running water.
 4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
 5. Place tip holder, clamping nut and rinsing tip in a suitable container into the ultra sonic device or instrument bath filled with detergent/disinfectant.
 6. Please note that, due to their internal office, the chemical disinfection of sonic and ultrasonic tips in an ultrasonic bath using a 2% solution always takes 10 minutes. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
 7. On completion of the immersion time, rinse tip holder, clamping nut and rinsing tip thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the tip holder and the inner geometry of the clamping nut with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage or the thread of the clamping nut.
 8. Dry tip holder, clamping nut and rinsing tip (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the tip holder and clamping nut.
 9. Visual examination to ensure that the tip holder, clamping nut and rinsing tip are clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.
 According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.
 Sterilisation in the autoclave:
 Make sure that the instruments is suitable for sonic tips classed as risk group critical B and the chosen method of sterilisation. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. In the sterilisation container: Place tip holder and clamping nut into a suitable sterilisation container (e.g. ref. 9952). For sterilisation, the container has to be sealed in suitable packaging. The rinsing tip must be sealed separately in a suitable sterilisation pouch.
 For this procedure, the following instructions have to be observed:
 Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
 -fractonated pre-vacuum (type B)
 -sterilisation temperature: 134°C
 -hold time: at least 5 minutes (full cycle)
 -drying time: at least 10 minutes
 In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standards DIN EN 13060. Make sure there is no exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.
 Universally valid notes:
 Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

8. Storage
 Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

9. Disposal
 The instruments have to be disposed of in: tight, non-breakable and non-pierceable containers (protection from contamination).

10. Safety and liability
 See also "Instructions for use and safety recommendations" in the current Komet Dental catalogue. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the appropriate use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to nonobservance of our recommendations for use or warnings as well as inadvertent misuse by the user.

11. Procedure in case of serious adverse events
 In case of serious adverse events, please proceed in line with the legal regulations in your country.

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© Gebr. Brasseler GmbH & Co. KG