



Optimal irrigation. Even better in blue.

BlueTip, the new sound-activated rinsing tip by Komet, makes the irrigation of root canals simple, fast and safe.

Lemgo, May 18, 2026 – In endodontics, cleaning and disinfecting the root canal system with as little residue as possible is a crucial factor for the success of the treatment. Since bacteria and the tissue that forms their breeding ground can be reduced very effectively by the sound-activated rinsing tip. Komet makes it easy, quick, and safe for users of BlueTip to achieve better cleaning results.

Optimized contrast for easier cleaning

BlueTip, the flexible rinsing tip made of polyamide, is more than just another addition to the range of perfectly coordinated high-quality instruments provided by the medical device manufacturer located in Lemgo, Germany. Above all, the blue color of the instrument significantly increases visibility in the canal as it contrasts clearly with the dentin.

The blue color of the BlueTip facilitates the insertion of the tip into the canal without bending it. During the rinsing process, it is easier to see whether the tip is oscillating freely or touching the canal walls. In addition, adhering debris or remnants of a medicated insert are more clearly visible and can be removed if necessary. This reduces blockages and makes it easier to remove residues.

This is all the more important as a root canal that is as residue-free as possible supports successful obturation and decisively contributes to the long-term success of the treatment. A positive side effect is that the BlueTip is easy to recognize on the tray and is always ready to hand.

Save time with "sonic speed"

BlueTip can be conveniently inserted into a standard sonic handpiece and activates the rinsing solution in the root canal. The solution activated in this manner requires a shorter reaction time to be effective, thus saving dentists valuable time. At the same time, the cleaning performance is enhanced. Bacteria, tissue and the smear layer are effectively reduced.

Follows the course of the canal precisely - flexible and delicate

BlueTip is made of polyamide to ensure optimum protection of natural tooth substance. Made of flexible synthetic material, the delicate tip follows the course of the canal precisely. This makes it noticeably easier for dentists to minimize the risk of undesirable material removal. In this way, BlueTip helps you to preserve the desired state of the canal anatomy.

| |
|--------------|
| Date |
| May 18, 2026 |

| |
|----------------------------------|
| Reprint restrictions |
| reprint free, please submit copy |

| |
|---|
| Contact Komet Dental |
| <p>Komet Dental Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 · 32657 Lemgo Postfach 160 · 32631 Lemgo Germany</p> <p>+49 (0) 5261 701-0 info@kometdental.de www.kometdental.com</p> <p>Limited Partnership, based in Lemgo Registration court Lemgo HRA 1156</p> <p>General partner: DENTAL Brasseler GmbH, based in Lemgo, Registration court Lemgo HRB 1157</p> <p>Managing directors: Stephan Köhler, Alexander Bartel, Markus Trachternach</p> |

| |
|--|
| Contact Press |
| <p>Julia Lempe Global Marketing Manager</p> <p>+49 (0) 5261 701-1890 jlempe@kometdental.com</p> |

Press Release



According to users, the fine polyamide tip fits into almost any canal. Also convenient: there is only one size to suit any application and, thanks to sterile single-use packaging units, there is no need for documentation. This makes integration into everyday practice and beyond extremely easy.

BlueTip complements Komet's endodontic portfolio, so that the complete range of instruments for root canal treatment is now perfectly coordinated and available in the high quality you have come to expect from Komet.

The following Images are available:

Image 1: BlueTip with sonic handpiece - close



Source: Komet Dental

Image 2: BlueTip with sonic handpiece



Source: Komet Dental

| |
|--------------|
| Date |
| May 18, 2026 |

| |
|----------------------------------|
| Reprint restrictions |
| reprint free, please submit copy |

| |
|---|
| Contact Komet Dental |
| <p>Komet Dental Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 · 32657 Lemgo Postfach 160 · 32631 Lemgo Germany</p> <p>+49 (0) 5261 701-0 info@kometdental.de www.kometdental.com</p> <p>Limited Partnership, based in Lemgo Registration court Lemgo HRA 1156</p> <p>General partner: DENTAL Brasseler GmbH, based in Lemgo, Registration court Lemgo HRB 1157</p> <p>Managing directors: Stephan Köhler, Alexander Bartel, Markus Trachternach</p> |

| |
|--|
| Contact Press |
| <p>Julia Lempe Global Marketing Manager</p> <p>+49 (0) 5261 701-1890 jlempe@kometdental.com</p> |



Image 3: BlueTip



Source: Komet Dental

About Komet Dental

Since 1923, Gebr. Brasseler GmbH & Co KG, headquartered in Lemgo, Germany, has been offering one of the world's largest ranges of rotary dental instruments and systems under the Komet Dental brand, with over 25,000 items. These are available through subsidiaries and partners in around 100 countries worldwide.

The comprehensive range includes sophisticated standard instruments and innovative solutions that set industry-wide standards in dental practices and laboratories. In addition to the classic "rose bur", the broad portfolio includes many other instruments such as crown cutters, composite removers and polishers. This is supplemented by products and systems for endodontic treatments. Services such as the Komet Academy with training courses on the use of the products complete the range.

About Brasseler

As a family business, Brasseler has been committed to improving the quality of life of people around the world with its dental and surgical products for over 100 years. What began in 1923 with a visionary idea has become an international success story: Having developed into a globally active MedTech company based in East Westphalia-Lippe, Brasseler is now one of the most important employers in the region. Brasseler has been family-owned since its foundation, and the shared values of the family business are characterized by long-term thinking and a focus on people. The company employs around 1,500 staff worldwide, 1,200 of whom work at the company head office in Lemgo.

| |
|--------------|
| Date |
| May 18, 2026 |

| |
|----------------------------------|
| Reprint restrictions |
| reprint free, please submit copy |

| |
|--|
| Contact Komet Dental |
| Komet Dental Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 · 32657 Lemgo Postfach 160 · 32631 Lemgo Germany +49 (0) 5261 701-0 info@kometdental.de www.kometdental.com Limited Partnership, based in Lemgo Registration court Lemgo HRA 1156 General partner: DENTAL Brasseler GmbH, based in Lemgo, Registration court Lemgo HRB 1157 Managing directors: Stephan Köhler, Alexander Bartel, Markus Trachternach |

| |
|---|
| Contact Press |
| Julia Lempe Global Marketing Manager +49 (0) 5261 701-1890 jlempe@kometdental.com |